

System II+ Cassettes



Optimised print surface improves scannability

CellPath's System II+ Cassettes come with a separate white lid that is easily attached via the use of a hinge. The robust locking mechanism on the cassettes close with a positive click and ensures the security of the specimen during processing.

There are multiple versions available; slotted, biopsy, micro biopsy and round holes, all of which have been validated for use with the Leica IPC, Sakura AutoWrite, EpreDia™ MicroWriter, Primera, and Vogel Printers. The laser versions can be used with Fa-Tech® cassette printers.



- A Slotted** – Incorporating 1.1mm wide slots which increases reagent flow, facilitating improved consistency and quality of processed specimens.
- B Biopsy** – Incorporating 0.9 by 0.9mm apertures, offering improved specimen security for small biopsy and tissue fragments.
- C Micro Biopsy** – Incorporating an inner chamber with 0.38 by 0.38mm apertures. The chamber's rounded corners allow for easier handling of samples using forceps.
- D Round** – Incorporating round 2.4mm diameter holes, facilitating increased reagent flow in a traditional cassette.

KEY BENEFITS

- 1 Verified for use with multiple cassette printers
- 2 Precisely engineered printing surface optimises print quality and scannability / read rates
- 3 Robust and secure closure ensures specimen security
- 4 Designed to improve reagent flow
- 5 Large colour range. Bespoke colours also available
- 6 Writing surface suitable for marking with a CellMark pen



Optimal quality printing of text and data matrix barcodes.

Each cassette is manufactured with a precisely engineered surface that optimises the print quality of text and linear / 2D data matrix barcodes. This reduces potential human errors such as misreading poorly printed or written text. It also significantly improves the scan rates, increasing the laboratory's efficiency, whilst reducing the need for human data entry and potential frustration.


We offer the cassettes in 15 colours, all a shade that provides greater contrast to printed text and barcodes which again improves the scannability and read rates. Bespoke colours are also available on request.

Tissue Processing Cassettes

The slotted, biopsy and round holed versions possess a drain hole to the rear of the cassette, allowing residual reagent trapped in the cassette when upright in a processing rack to better drain and reduce carry over and spoiling of reagents.

The micro biopsy version incorporates eight sink compartments which enable reagents to flow more freely into the cassette and air to be expelled, which might ordinarily have been trapped by the small apertures.

This further improves the reagent flow.



All versions are supplied loose, threaded or hopped giving customers a wide choice of options depending on their requirements. The threaded and hopped options allow multiple cassettes to be easily loaded into a cassette printer at once, saving considerable laboratory staff time.



The System II+ Cassettes are available either loose, in hopper tubes or threaded.

The System II+ Cassettes are available in a full range of colours



CellPath Ltd
80 Mochdre Enterprise Park,
Newtown, Powys, SY16 4LE, UK
T: +44 (0)1686 611333
E: sales@cellpath.com
cellpath.com



CARING FOR OUR
PLANET

Helping our customers meet their environmental commitments
















CellPath are committed to making the world a greener place and keeping our ecological footprint as small as possible.

As such, cassette threads are now manufactured from

specialist cardboard made from pulp that is recyclable and fully biodegradable.
















Our hopper tubes are made from rPet (*recycled polyethylene terephthalate*) which is fully recyclable.


Inkjet/Thermal Printer Compatible System II+ Slotted Cassettes

Loose		Hoppered		Threaded		Colour
Item Code	Quantity	Item Code	Quantity	Item Code	Quantity	
EAM-0502-03A	1000	EAM-0502-05H	3000	EAM-0502-72B	2000	White 
EAM-0504-03A	1000	EAM-0504-05H	3000	EAM-0504-72B	2000	Green 
EAM-0505-03A	1000	EAM-0505-05H	3000	EAM-0505-72B	2000	Grey 
EAM-0506-03A	1000	EAM-0506-05H	3000	EAM-0506-72B	2000	Blue 
EAM-0507-03A	1000	EAM-0507-05H	3000	EAM-0507-72B	2000	Orange 
EAM-0508-03A	1000	EAM-0508-05H	3000	EAM-0508-72B	2000	Pink 
EAM-0509-03A	1000	EAM-0509-05H	3000	EAM-0509-72B	2000	Yellow 
EAM-0510-03A	1000	EAM-0510-05H	3000	EAM-0510-72B	2000	Red 
EAM-0511-03A	1000	EAM-0511-05H	3000	EAM-0511-72B	2000	Peach 
EAM-0512-03A	1000	EAM-0512-05H	3000	EAM-0512-72B	2000	Tan 
EAM-0515-03A	1000	EAM-0515-05H	3000	EAM-0515-72B	2000	Brown 
EAM-0516-03A	1000	EAM-0516-05H	3000	EAM-0516-72B	2000	Dark Blue 
EAM-0518-03A	1000	EAM-0518-05H	3000	EAM-0518-72B	2000	Dark Green 
EAM-0519-03A	1000	EAM-0519-05H	3000	EAM-0519-72B	2000	Lavender 
EAM-0520-03A	1000	EAM-0520-05H	3000	EAM-0520-72B	2000	Turquoise 


















Inkjet/Thermal Printer Compatible System II+ Biopsy Cassettes

Loose		Hoppered		Threaded		Colour
Item Code	Quantity	Item Code	Quantity	Item Code	Quantity	
EAM-0602-03A	1000	EAM-0602-05H	3000	EAM-0602-72B	2000	White 
EAM-0604-03A	1000	EAM-0604-05H	3000	EAM-0604-72B	2000	Green 
EAM-0605-03A	1000	EAM-0605-05H	3000	EAM-0605-72B	2000	Grey 
EAM-0606-03A	1000	EAM-0606-05H	3000	EAM-0606-72B	2000	Blue 
EAM-0607-03A	1000	EAM-0607-05H	3000	EAM-0607-72B	2000	Orange 
EAM-0608-03A	1000	EAM-0608-05H	3000	EAM-0608-72B	2000	Pink 
EAM-0609-03A	1000	EAM-0609-05H	3000	EAM-0609-72B	2000	Yellow 
EAM-0610-03A	1000	EAM-0610-05H	3000	EAM-0610-72B	2000	Red 
EAM-0611-03A	1000	EAM-0611-05H	3000	EAM-0611-72B	2000	Peach 
EAM-0612-03A	1000	EAM-0612-05H	3000	EAM-0612-72B	2000	Tan 
EAM-0615-03A	1000	EAM-0615-05H	3000	EAM-0615-72B	2000	Brown 
EAM-0616-03A	1000	EAM-0616-05H	3000	EAM-0616-72B	2000	Dark Blue 
EAM-0618-03A	1000	EAM-0618-05H	3000	EAM-0618-72B	2000	Dark Green 
EAM-0619-03A	1000	EAM-0619-05H	3000	EAM-0619-72B	2000	Lavender 
EAM-0620-03A	1000	EAM-0620-05H	3000	EAM-0620-72B	2000	Turquoise 



















Inkjet/Thermal Printer Compatible System II+ Micro Biopsy Cassettes

Loose		Hoppered		Threaded		Colour
Item Code	Quantity	Item Code	Quantity	Item Code	Quantity	
EAM-0702-03A	1000	EAM-0702-05H	3000	EAM-0702-72B	2000	White 
EAM-0704-03A	1000	EAM-0704-05H	3000	EAM-0704-72B	2000	Green 
EAM-0705-03A	1000	EAM-0705-05H	3000	EAM-0705-72B	2000	Grey 
EAM-0706-03A	1000	EAM-0706-05H	3000	EAM-0706-72B	2000	Blue 
EAM-0707-03A	1000	EAM-0707-05H	3000	EAM-0707-72B	2000	Orange 
EAM-0708-03A	1000	EAM-0708-05H	3000	EAM-0708-72B	2000	Pink 
EAM-0709-03A	1000	EAM-0709-05H	3000	EAM-0709-72B	2000	Yellow 
EAM-0710-03A	1000	EAM-0710-05H	3000	EAM-0710-72B	2000	Red 
EAM-0711-03A	1000	EAM-0711-05H	3000	EAM-0711-72B	2000	Peach 
EAM-0712-03A	1000	EAM-0712-05H	3000	EAM-0712-72B	2000	Tan 
EAM-0715-03A	1000	EAM-0715-05H	3000	EAM-0715-72B	2000	Brown 
EAM-0716-03A	1000	EAM-0716-05H	3000	EAM-0716-72B	2000	Dark Blue 
EAM-0718-03A	1000	EAM-0718-05H	3000	EAM-0718-72B	2000	Dark Green 
EAM-0719-03A	1000	EAM-0719-05H	3000	EAM-0719-72B	2000	Lavender 
EAM-0720-03A	1000	EAM-0720-05H	3000	EAM-0720-72B	2000	Turquoise 



















Inkjet/Thermal Printer Compatible System II+ Round Cassettes

Loose		Hoppered		Threaded		Colour
Item Code	Quantity	Item Code	Quantity	Item Code	Quantity	
EAM-0802-03A	1000	EAM-0802-05H	3000	EAM-0802-72B	2000	White 
EAM-0804-03A	1000	EAM-0804-05H	3000	EAM-0804-72B	2000	Green 
EAM-0805-03A	1000	EAM-0805-05H	3000	EAM-0805-72B	2000	Grey 
EAM-0806-03A	1000	EAM-0806-05H	3000	EAM-0806-72B	2000	Blue 
EAM-0807-03A	1000	EAM-0807-05H	3000	EAM-0807-72B	2000	Orange 
EAM-0808-03A	1000	EAM-0808-05H	3000	EAM-0808-72B	2000	Pink 
EAM-0809-03A	1000	EAM-0809-05H	3000	EAM-0809-72B	2000	Yellow 
EAM-0810-03A	1000	EAM-0810-05H	3000	EAM-0810-72B	2000	Red 
EAM-0811-03A	1000	EAM-0811-05H	3000	EAM-0811-72B	2000	Peach 
EAM-0812-03A	1000	EAM-0812-05H	3000	EAM-0812-72B	2000	Tan 
EAM-0815-03A	1000	EAM-0815-05H	3000	EAM-0815-72B	2000	Brown 
EAM-0816-03A	1000	EAM-0816-05H	3000	EAM-0816-72B	2000	Dark Blue 
EAM-0818-03A	1000	EAM-0818-05H	3000	EAM-0818-72B	2000	Dark Green 
EAM-0819-03A	1000	EAM-0819-05H	3000	EAM-0819-72B	2000	Lavender 
EAM-0820-03A	1000	EAM-0820-05H	3000	EAM-0820-72B	2000	Turquoise 


















Laser Printer Compatible System II+ Slotted Cassettes

Threaded		
Item Code	Quantity	Colour
EAM-1302-72B	2000	White 
EAM-1304-72B	2000	Green 
EAM-1305-72B	2000	Grey 
EAM-1306-72B	2000	Blue 
EAM-1307-72B	2000	Orange 
EAM-1308-72B	2000	Pink 
EAM-1309-72B	2000	Yellow 
EAM-1310-72B	2000	Red 
EAM-1311-72B	2000	Peach 
EAM-1312-72B	2000	Tan 
EAM-1315-72B	2000	Brown 
EAM-1316-72B	2000	Dark Blue 
EAM-1318-72B	2000	Dark Green 
EAM-1319-72B	2000	Lavender 
EAM-1320-72B	2000	Turquoise 


















Laser Printer Compatible System II+ Biopsy Cassettes

Threaded		
Item Code	Quantity	Colour
EAM-1402-72B	2000	White 
EAM-1404-72B	2000	Green 
EAM-1405-72B	2000	Grey 
EAM-1406-72B	2000	Blue 
EAM-1407-72B	2000	Orange 
EAM-1408-72B	2000	Pink 
EAM-1409-72B	2000	Yellow 
EAM-1410-72B	2000	Red 
EAM-1411-72B	2000	Peach 
EAM-1412-72B	2000	Tan 
EAM-1415-72B	2000	Brown 
EAM-1416-72B	2000	Dark Blue 
EAM-1418-72B	2000	Dark Green 
EAM-1419-72B	2000	Lavender 
EAM-1420-72B	2000	Turquoise 

















Tissue Processing Cassettes

Laser Printer Compatible System II+ Micro Biopsy Cassettes

Threaded		
Item Code	Quantity	Colour
EAM-1502-72B	2000	White 
EAM-1504-72B	2000	Green 
EAM-1505-72B	2000	Grey 
EAM-1506-72B	2000	Blue 
EAM-1507-72B	2000	Orange 
EAM-1508-72B	2000	Pink 
EAM-1509-72B	2000	Yellow 
EAM-1510-72B	2000	Red 
EAM-1511-72B	2000	Peach 
EAM-1512-72B	2000	Tan 
EAM-1515-72B	2000	Brown 
EAM-1516-72B	2000	Dark Blue 
EAM-1518-72B	2000	Dark Green 
EAM-1519-72B	2000	Lavender 
EAM-1520-72B	2000	Turquoise 



Laser Printer Compatible System II+ Round Cassettes

Threaded		
Item Code	Quantity	Colour
EAM-1602-72B	2000	White 
EAM-1604-72B	2000	Green 
EAM-1605-72B	2000	Grey 
EAM-1606-72B	2000	Blue 
EAM-1607-72B	2000	Orange 
EAM-1608-72B	2000	Pink 
EAM-1609-72B	2000	Yellow 
EAM-1610-72B	2000	Red 
EAM-1611-72B	2000	Peach 
EAM-1612-72B	2000	Tan 
EAM-1615-72B	2000	Brown 
EAM-1616-72B	2000	Dark Blue 
EAM-1618-72B	2000	Dark Green 
EAM-1619-72B	2000	Lavender 
EAM-1620-72B	2000	Turquoise 